



**ROTAID™**  
AED CABINETS



**ROTAID LIST**  
**AED Overview**



**Serial number**

# Rotaid Cabinet

The serial number of the Rotaid AED cabinet can be found on your invoice and on the location shown here:



# SOLID PLUS & SOLID PLUS HEAT LED AED PLACEMENT ASSISTANT





Place the defibrillator (AED) on the platform.  
Use the elastic band to secure the defibrillator.

<b>Cardiac Science</b> Powerheart G5	<b>Nihon Kohden</b> Cardiolife AED-2150 Cardiolife AED-3100	<b>Radian</b> HR-501 HR-503 HR-503A
<b>CU Medical Systems</b> I-Pad CU-SP1 I-Pad NFI200	<b>Philips</b> Heartstart FR3 Heartstart FRx Heartstart OnSite/HS1	<b>Schiller</b> Fred Easy Fred Easy Life Fred Easyport Fred PA-1
<b>Defibtech</b> Lifeline Pro Lifeline View	<b>Physio control</b> CR2 Lifepak 1000 Lifepak CR Plus Lifepak Express	<b>Zoll</b> AED 3 AED Plus
<b>HeartSine</b> Samaritan Pad 350P Samaritan Pad 360P Samaritan Pad 450P	<b>Primedic</b> Heartsave AED Heartsave AS Heartsave One Heartsave Pad	
<b>Mediana</b> A10 A15	<b>Progetti</b> Rescue Sam Rescue Sam Pro	
<b>Mindray</b> BeneHeart D1 BeneHeart D1 Pro		

# SOLID PLUS & SOLID PLUS HEAT LED AED PLACEMENT ASSISTANT

Use 'both bolts' (2x DIN screw M6) to hang the AED.  
If needed, the AED can also be placed on a 'special bracket'.



<b>Cardiac Science</b> Powerheart G3 Plus Powerheart G3 Pro Powerheart G5	<b>Mindray</b> BeneHeart D1 BeneHeart D1 Pro	<b>Progetti</b> Rescue Sam Rescue Sam Pro
<b>Cardiaid</b> CT0207RF CT0207RS	<b>Nihon Kohden</b> Cardiolife AED-2150 Cardiolife AED-3100	<b>Radian</b> HR-501 HR-503 HR-503A
<b>CU Medical Systems</b> I-Pad CU-SP1 I-Pad NF1200	<b>Philips</b> Heartstart FR3 Heartstart FRx Heartstart OnSite/HS1	<b>Schiller</b> Fred Easy Fred Easy Life Fred Easyport Fred PA-1
<b>Defibtech</b> Lifeline Lifeline Pro Lifeline View	<b>Physio control</b> CR2 Lifepak 1000 Lifepak CR Plus Lifepak Express	<b>Zoll</b> AED 3 AED Plus
<b>HeartSine</b> Samaritan Pad 350P Samaritan Pad 360P Samaritan Pad 450P	<b>Primedic</b> Heartsave AED Heartsave AS Heartsave One Heartsave Pad	
<b>Mediana</b> A10 A15		

# ROTAID AED LIST

## AED List

\* Click on your AED and go directly to the correct page:



### Cardiac Science

Powerheart G3 Plus  
Powerheart G3 Pro  
Powerheart G5



### Philips

Heartstart FR3  
Heartstart FRx  
Heartstart OnSite/HS1



### Cardiaid

CT0207RF  
CT0207RS



### Physio control

CR2  
Lifepak 1000  
Lifepak CR Plus  
Lifepak Express



Medical Systems, Inc.

### CU Medical Systems

I-Pad CU-SP1  
I-Pad NF1200



### Defibtech

Lifeline  
Lifeline Pro  
Lifeline View



### Primedic

Heartsave AED  
Heartsave AS  
Heartsave One  
Heartsave Pad



### HeartSine

Samaritan Pad 350P  
Samaritan Pad 360P  
Samaritan Pad 450P



### Progetti

Rescue Sam  
Rescue Sam Pro



### Mediana

A10  
A15



### Schiller

Fred Easy  
Fred Easy Life  
Fred Easyport  
Fred PA-1



### Mindray

BeneHeart D1  
BeneHeart D1 Pro



### Telefunken

Telefunken AED



### Nihon Kohden

Cardiolife AED-2150  
Cardiolife AED-3100



### Zoll (install without bag)

AED 3  
AED Plus



Powerheart G3 (Elite/Plus)



Powerheart G5





**CardiAid<sup>®</sup>**

CT0207R (F/S)





Medical Systems, Inc.

## I-Pad (CU-SP1/NF1200)







## Lifeline (Pro)



## Lifeline View





## Samaritan PAD (350P/360P/450P)





A10



A15





**mindray**

## Beneheart D1 (Pro)





CardioLife AED-2100



CardioLife AED-3100





# PHILIPS

## Heartstart (HS1/FRx)



HS1

FRx





CR2

Lifepak (1000/CR Plus/Express)





### Heartsave (Aed/AS)



### Heartsave (One/Pad)







## Rescue Sam (PRO)





**SCHILLER**

The Art of Diagnostics

Fred Easy (Life/-port)





**TELEFUNKEN**

Telefunken





**ZOLL**

**AED 3** (install without bag)

**AED Plus** (install without bag)

